

OmniOx (HFT 500)



Jan. 2013. Prepared by MEKICS



1. Introduce of OmniOx (HFT500)
2. Specification (HFT, CPAP+)
3. The World Best Features



1. Introduce

1. What is **High Flow Therapy's** Definition?

- ✓ Flows that exceed patient demands at various minute volumes

2. How does HFT impact breathing?

- CO2 Ventilation
 - ✓ Washout of the nasopharyngeal dead space
- Efficient Oxygenation
- Work of Breathing
 - ✓ Reduction in inspiratory resistance associated with gas flow through the nasopharynx
- Energy Cost of Gas Conditioning
 - ✓ Improvement in respiratory mechanical parameters associated with gas temperature and state of humidification
 - ✓ Reduction in metabolic work associated with gas conditioning
- Provision of mild distending pressure




1. Introduce

3. Must reduce the dead space of Upper Airway (the nasal cavity)


- ✓ Circulation and Maintain of fresh air (minimized dead space)
- ✓ Appropriated PEEP (airway security)
- ✓ Minimized resistance of exhalation (efficiency of CO₂ deflation and optimized positive pressure)

→ The best solution is
"JET FLOW EFFECT"



4. CO₂ Ventilation

- ✓ Washout of nasopharyngeal dead space
- ✓ Improved fractions of alveolar gases with respect to carbon dioxide
- ✓ Low flow nasal cannula therapy is only thought to facilitate oxygenation
- ✓ HFT impacts CO₂ elimination





Flushing of dead space in the Nasopharyngeal cavity helps Enhance alveolar ventilation



1. Introduce

5. Efficient Oxygenation

- ✓ High inspired oxygen fractions by eliminating room air entrainment during inspiration
- ✓ Improve alveolar oxygen fractions beyond mask therapy based on the equation for alveolar ventilation
- ✓ Patients can often maintain better oxygenation or require a lower FiO_2 compared to conventional mask or cannula therapies



6. Work of Breathing(WOB)

- ✓ HFT provides enough flow to match or exceed a patient's inspiratory flow
- ✓ HFT most likely minimizes the inspiratory resistance associated with the nasopharynx
- ✓ **Adequate warming and humidification** of the conducting airways by delivery of warm, humid gas is associated with improved conductance and pulmonary compliance compared to dry, cooler gas
- ✓ Delivery of breathing gases at body temperature and saturation **promotes an ideal respiratory mechanical response**



7. Energy Cost of Gas Conditioning

- ✓ The nasal air passages expend energy to warm inspiratory air from ambient to 37 °C and vaporize water to humidify the incoming air to 100% relative humidity
- ✓ Alleviated when gas is delivered at body temperature and saturated



8. ACPAP+ ?


- Unique technology for the first time by MEKICS in the World
 - ✓ New concept of positive pressure support for easy exhalation that support inspiration in inhalation of breath and interrupt the flow quickly in exhalation
 - ✓ Same effect but more comfortable to inhale and exhale both together so it is basically big different with CPAP.
- Clinical effect of comfortable breathing
 - ✓ Most of case of CPAP patient needs to have adaptation period due to uncomfortable
 - ✓ Increase of WOB, and disturb of stable sleeping
 - ✓ Clinical effect beside the ventilation only.



1. Introduce

ACPAP+® (Adaptive CPAP & PSV)


- Positive Pressure Support, Lower expiration resistance as support of inspiration
- Maintain of airway and comfortable breathing



- ✓ Easy to breathing?
= Possible to support of inspiration and low of expiration resistance
- ✓ Support of Inspiration: PSV Trg
- ✓ Expiration Resistance : Ex_Sen, ExTime



HFT500 Target Market



2. Specification

● Functions

➤ General

- ✓ 4.3" Color TFT, Touch Screen, Knob
- ✓ Mode : HF/CPAP, PS/CPAP
- ✓ O2 Mixer
- ✓ **FiO2, Respiration, SpO2 Monitoring**
- ✓ Built in Humidifier

➤ HF/CPAP : Flow(Lpm), O2(%)

➤ PS/CPAP : PSV, PAP, Trigger, ExTime, Ex_Sense, O2(%), Wake-up

● Performance



➤ Pressure : Up to 100cmH2O

➤ Flow : Up to 160 Lpm / 60 Lpm (HFT)

➤ O2 : 21 ~ 100%



3. The World Best Features (HF/CPAP)





	MEKICS	Competitors
Flow /Pressure Drive Power	10~ 60 lpm Up to 100cmH2O	10 ~ 40 lpm(60lpm) Up to 40cmH2O
Effect	Available Jet Flow Low Expiratory Resistance	Not Jet-Flow High Expiratory Resistance

3. The World Best Features (PS/CPAP)


ACPAP+®

PSV at Ins. Phase

Quick Pause in Exp. Phase



<MEKICS>



<Competitors>



References

1. **Finer NN, Bates R, Tomat P.** Low flow oxygen delivery via nasal cannula to neonates. *Pediatr Pulmonol* 1996;21(1):48-51.
2. **Sreenan C, Lemke RP, Hudson-Mason A, Osiovich H.** High-flow nasal cannulae in the management of apnea of prematurity: a comparison with conventional nasal continuous positive airway pressure. *Pediatrics* 2001;107(5):1081-1083.
3. **Waugh JB, Granger WM.** An evaluation of 2 new devices for nasal high-flow gas therapy. *Respir Care* 2004;49(8):902-906.
4. **Negus VE.** Humidification of the air passages. *Thorax* 1952;7(2):148-151.
5. **Williams R, Rankin N, Smith T, Galler D, Seakins P.** Relationship between the humidity and temperature of inspired gas and the function of the airway mucosa. *Crit Care Med* 1996;24(11):1920-1929.
6. **Kelly MG, McGarvey LP, Heaney LG, Elborn JS.** Nasal septal perforation and oxygen cannulae. *Hosp Med* 2001;62(4):248.
7. **Robertson NJ, McCarthy LS, Hamilton PA, Moss AL.** Nasal deformities resulting from flow driver continuous positive airway pressure. *Arch Dis Child Fetal Neonatal Ed* 1996;75(3):F209-212.
8. **Loftus BC, Ahn J, Haddad J, Jr.** Neonatal nasal deformities secondary to nasal continuous positive airway pressure. *Laryngoscope* 1994;104(8 Pt 1):1019-1022.
9. **Kopelman AE, Holbert D.** Use of oxygen cannulas in extremely low birthweight infants is associated with mucosal trauma and bleeding, and possibly with coagulase-negative staphylococcal sepsis. *J Perinatol* 2003;23(2):94-97.
10. **Campbell EJ, Baker MD, Crites-Silver P.** Subjective effects of humidification of oxygen for delivery by nasal cannula. A prospective study. *Chest* 1988;93(2):289-293.
11. **Rankin N.** What is optimum humidity? *Respir Care Clin N Am* 1998;4(2):321-328.
12. **Williams RB.** The effects of excessive humidity. *Respir Care Clin N Am* 1998;4(2):215-228.
13. **Woodhead DD, Lambert DK, Clark JM, Christensen RD.** Comparing two methods of delivering high-flow gas therapy by nasal cannula following endotracheal extubation: a prospective, randomized, masked, crossover trial. *J Perinatol* 2006;26(8):481-485.
14. **Dewan NA, Bell CW.** Effect of low flow and high flow oxygen delivery on exercise tolerance and sensation of dyspnea. A study comparing the transtracheal catheter and nasal prongs. *Chest* 1994;105(4):1061-1065.



References

15. **Shepard JW, Jr., Burger CD.** Nasal and oral flow-volume loops in normal subjects and patients with obstructive sleep apnea. *Am Rev Respir Dis* 1990;142(6 Pt 1):1288-1293.
16. **Greenspan JS, Wolfson MR, Shaffer TH.** Airway responsiveness to low inspired gas temperature in preterm neonates. *J Pediatr* 1991;118(3):443-445.
17. **Fontanari P, Burnet H, Zattara-Hartmann MC, Jammes Y.** Changes in airway resistance induced by nasal inhalation of cold dry, dry, or moist air in normal individuals. *J Appl Physiol* 1996;81(4):1739-1743.
18. **Fontanari P, Zattara-Hartmann MC, Burnet H, Jammes Y.** Nasal eupnoeic inhalation of cold, dry air increases airway resistance in asthmatic patients. *Eur Respir J* 1997;10(10):2250-2254.
19. **Proctor DF.** Physiology of the upper airway. In: Visher MB, Hastings AB, Pappenheimer JR, Rahn H, editors. *Handbook of Physiology-Respiration 1*. Baltimore: Williams & Wilkins; 1985. p 309-345.
20. **Mlynski G.** Physiology and pathophysiology of nasal breathing. In: Behrbohm H, Tardy T, editors. *Essentials of Septorhinoplasty: Philosophy-Approaches-Techniques*. Stuttgart, NY: Thieme Medical Publishers; 2004. p 75-87.
21. **Locke RG, Wolfson MR, Shaffer TH, Rubenstein SD, Greenspan JS.** Inadvertent administration of positive end-distending pressure during nasal cannula flow. *Pediatrics* 1993;91(1):135-138.
22. **Kahn DJ, Courtney SE, Steele AM, Habib RH.** Unpredictability of Delivered Bubble Nasal Continuous Positive Airway Pressure Role of Bias Flow Magnitude and Nares-Prong Air Leaks. *Pediatr Res* 2007.
23. **Kubicka ZJ, Limauro J, Darnall RA.** Heated, humidified high-flow nasal cannula therapy: yet another way to deliver continuous positive airway pressure? *Pediatrics* 2008;121(1):82-88.
24. **Wilkinson DJ, Andersen CC, Smith K, Holberton J.** Pharyngeal pressure with high-flow nasal cannulae in premature infants. *J Perinatol* 2008;28(1):42-47.
25. **Saslow JG, Aghai ZH, Nakhla TA, Hart JJ, Lawrysh R, Stahl GE, Pyon KH.** Work of breathing using high-flow nasal cannula in preterm infants. *J Perinatol* 2006;26(8):476-480.
26. **Pyon KH, Aghai ZH, Nakhla TA, Stahl GE, Saslow JG.** High Flow Nasal Cannula in Preterm Infants: Effects of High Flow Rates on Work of Breathing. *PAS* 2008;63:3763.3713.
27. **Spence KL, Murphy D, Kilian C, McGonigle R, Kilani RA.** High-flow nasal cannula as a device to provide continuous positive airway pressure in infants. *J Perinatol* 2007;27(12):772-775.



END

Thanks for your attention

