OPERATION VIDEO/ANLEITUNGSVIDEO/VIDEO TUTORIAL/TUTORIEL VIDEO/VIDEO TUTORIAL

Instruction for SARS-CoV-2
Antigen Test Kit (Colloidal Gold)

REF: 52104115

- For self-testing
- For use by people between 18 and 75 years of age. For child and adolescent ages 18 years or younger and for people over 75 years of age, use only under the supervision of authorized persons.
- PRECAUTIONS**
- Please read all the information in this instruction before performing the test.**
- This kit is only used for in vitro diagnosis.**
- The test card is for one-time use only, and cannot be reused or used by more than one person at a time.**
- Do not interchange components from different batches.**
- Do not use expired products. If the product packaging bag is damaged, poorly sealed or the sample diluent leaks, the product should not be used.**
- Avoid high temperatures (above 30°C) during the test. Cryogenically stored test cards, samples and sample diluents must be brought back to room temperature (18-28°C) prior to opening. Please do not perform the test in the following situations or scenarios: outdoors under the sunlight, temperatures above 30°C, in a moving vehicle or in any other place where stability cannot be guaranteed.**
- Do not touch the reaction area of the test paper and keep the swab clean.**
- For substances that contain or are suspected to contain sources of infection, appropriate biosafety assurance procedures should be followed. If the liquid splashes in your eyes, on your skin, rinse with plenty of water. If you feel unwell, go to a specialist for a checkup.**

SAMPLE REQUIREMENTS

- Nasal swab can be used for testing.

INTENDED USE & BENEFITS

Genrui SARS-CoV-2 antigen test kit (colloidal gold) is an immunochromatographic assay for rapid and qualitative test of N protein of Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) antigen detection in nasal swabs. This test is applicable to all individuals suspected of being infected by SARS-CoV-2. The test will be used to assist in the diagnosis of the coronavirus infectious disease (COVID-19) caused by SARS-CoV-2.

●The test kit is simple, safe, effective and intended for self-testing, which is suitable for individuals to use in non-laboratory settings like homes, offices, schools, sports stadiums, airports, etc.

TEST PRINCIPLE

During the test, the processed samples to be tested are added to the sample holes. When the sample contains SARS-CoV-2 antigen, it first combines with the colloidal gold-labelled anti-SARS-CoV-2 antigen antibody. Chromatography is then performed. When it binds to the anti-SARS-CoV-2 monoclonal antibody previously immobilized on another membrane, a purple-red band will appear in the test area (T). If there is no SARS-CoV-2 antigen in the sample, there will be no purple-red band in the test area (T). A purple-red band will appear in the quality control area (C) regardless of the presence of new coronavirus antigen in the sample, which is used as a criterion to determine whether there is sufficient sample or the chromatography is processed properly. It is also used as an internal control standard for the test kit.

MAIN COMPONENTS

- 1 test card, 1sample diluent (0.5 mL), 1 dripper, 1 swab, 1 biohazard sample bag and 1 instruction for use (IFU).
- The main components of the test card are foil pouch, desiccant, glass fiber membrane (colloidal gold labeled anti-novel coronavirus monoclonal antibody), nitrocellulose membrane (detection area coated with anti-novel coronavirus monoclonal antibody, quality control area coated with Goat anti Mouse IgG) and PVC plate.
- Sample diluent, the main component is phosphate buffer (PBS).
- Accessories required but not provided: Timer.

STORAGE CONDITIONS

- The test kit can be stored at 2-30°C and the expiration date is 18 months. Once the foil pouch is opened, the validity period is 1 hour at a temperature of 18-28°C and humidity less than 65%.
- The sample solution is valid for 1 month from the date of opening. The production date is shown on the outer box.

USE CONDITIONS

- Please make sure that the ambient temperature is 18-28°C when using.
- If the ambient humidity is below 65%, please use the product within one hour after opening the bag. And when the humidity is higher than 65%, please make sure to use the product immediately after opening it.

TRANSPORT CONDITIONS

Transport at 2-30°C.

QUALITY CONTROL

Each test card has a built-in control. The purple-red band at the control line can be considered as an internal positive program control. If the procedure is performed correctly, the control line will appear. If no control line appears, then the test is invalid and a new test should be performed.

LIMITATIONS

- Test results cannot be used as a basis for diagnosis. Comprehensive judgments should be made based on clinical symptoms, epidemiological conditions and further clinical data.
- The accuracy of the test depends on the sample collection process. Improper sample collection will affect the test results.

- A positive test result does not rule out concurrent infection with other pathogens. The negative result may be caused by:
 - a) Improper sample collection, improper sample transfer or handling so that the amount of virus in the sample is too low.
 - b) The level of SARS-CoV-2 antigen is lower than the limit of detection.
 - c) Variations in viral genes may lead to changes in antibody determinants.

- There may be other unlisted reasons that cause the detection to be abnormal.
- This product can not detect the SARS-CoV-2 antigen in the sample, but cannot determine the concentration of the antigen in the sample.

PERFORMANCE CHARACTERISTICS

Clinical performance

207 nasal swab samples, which include 107 confirmed as COVID-19 positive and 100 confirmed as COVID-19 negative by RT-PCR assay, were obtained for testing, and then compared the test results of Genrui SARS-CoV-2 Antigen Test Kit (Colloidal Gold) with that of RT-PCR. The results are shown below:

Positive Percent Agreement:

Ct value	Number of Sample	Number of true positive Rapid Test Samples	Number of false negative Rapid Test Samples	Sensitivity of SARS-CoV-2 Antigen Rapid Test (Ct)
≤30	80	80	0	100%(96-100)
≤32	91	91	0	100%(96-100)
≤34	101	100	1	99%(94-99)
≤36	107	105	2	98.13%(93-99)

Negative Percent Agreement:

Number of Sample	Number of true negative Rapid Test Samples	Number of false positive Rapid Test Samples	Specificity of SARS-CoV-2 Antigen Rapid Test (Ct)
100	100	0	100%(96-100)

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