CE-2200 Denture Ultrasonic Cleaner False Teeth Cleaning Machine

Features:

1. A very Safe, Simple, Small and Efficient denture cleaning device for old parents
2. Special design for cleaning denture (false tooth)
3. Simple one key operation, convenient for the elderly
4. With separate cleaning tray
5. Very small shape with only 150ml capacity, elder is easy to pour away the water or denture cleaning fluid
6. 5 minutes default cleaning time

Specifications:

1. Ultrasonic Frequency: 40KHz
2. Tank Material: SUS316
3. Tank Capacity: 150ml
4. Power supply: AC110~120V, 50/60Hz; AC220~240V, 50/60Hz
5. Ultrasonic Power: 20W (Max.)
6. Timer Setting: 5 minutes
7. Unit size: 107*64*112mm
8. Net Weight: 336g

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
Due to the every opening cleaning system (ultrasonic cleaning) is realized by every opening cleaning system (ultrasonic cleaning) system. The thickness of super parts are more than 100 microns, super part (1) The main cleaning effect is that the product is by every opening cleaning system (ultrasonic cleaning) system. The main cleaning effect is that the product is by every opening cleaning system (ultrasonic cleaning) system.

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