

Mini Fascia Gun Manual



④

⑤

⑥

⑦

⑧

Component Name:

Model: OM-MG03

① LED indicator light, ② power on/off/gear button, ③ Type-C charging port, ④ message gun body, ⑤ U-shaped head, ⑥ spherical head, ⑦ bullet shaped head, ⑧ flat head

Massage head suitable for specific areas:

- U-shaped head: suitable for the waist, back, shoulders, calves, and other areas.
- Spherical head: suitable for shoulders, waist, back, arms, hips, thighs, calves, and other large muscle groups.
- Bullet shaped head: suitable for the soles of the feet, palms, and other parts.
- Flat head: suitable for various muscle areas.

Operation Manual:

- Power on: Press and hold the power button for 2 seconds until the LED light comes on, indicating power on.
- Shutdown: Press and hold the power button for 2 seconds until the fascia gun stops working and the indicator light goes out, indicating that it has been shut down.
- Intensity gear switching: In the power on state, short press the power button to cycle between "1-6 gears", and the corresponding indicator light will remain on.
- After 15 minutes of continuous use, the device will automatically shut down.

Indicator light status description:

1.During work, the LED indicator light is constantly on. 2.During charging, the LED indicator light flashes. 3.When fully charged, the 6 LED indicator lights remain on.

Precautions:

1. The product has a certain amount of pre charge when it leaves the factory. It is recommended to fully charge it before use after unpacking. Do not use it while charging. After charging is complete, please unplug the charging cable.
2. Do not plug or unplug the charging cable with wet hands, otherwise it may cause electric shock and other accidents.
3. Patients with physical abnormalities are prohibited from using this product.
4. It is recommended to use it for no more than 30 minutes per day. Do not use continuously in the same area for more than 30 minutes.
- 5.After continuous use of this product, it is normal for the skin surface of the application area to show slight redness without any other discomfort.
6. Store this product in a dry and ventilated place, away from direct sunlight, high temperatures, and humidity. And prevent severe vibration, collision, rolling, and falling.
7. It is prohibited to use it in gas stations, high dust workshops, and gas stations with flammable and explosive hazards. Before cleaning this product, please make sure to disconnect the power supply. To clean the surface of the product, you can wipe it with a damp cloth instead of rinsing it directly with water.
9. This product contains lithium-ion batteries. Please do not disassemble, modify, immerse the product in water or put it into fire at will.
- 10.When disposing of the product, the power should be turned off, the lithium battery should be removed, and it should be safely disposed of according to the government regulations at that time.

Warning: This device can be used by children aged 8 and above, as well as individuals with weakened physical, sensory, or mental abilities or lacking experience and knowledge, provided that they have received supervision or guidance on safe use of the device and are aware of related hazards. Children are not allowed to play with this device. Cleaning and user maintenance must not be carried out by children without supervision.