



MICRO INFUSION KIT

User Manual

PLEASE READ

Before using the device, please ensure that you have carefully read and fully understood the user manual.

To use the device properly, make sure the instructions are followed carefully.

If you have any questions while using the device, please reach out to our company for assistance.

Safety Warnings:

- When using this product, strictly follow the instructions in the manual.
- Do not use the device in any method other than as the instructions in this user manual.
- Strictly adhere to the instructions for the use, sanitisation, and storage of the device and accessories.
- Use only the accessories provided by our company.

This product is classified as a cosmetic device for home use and is not intended to diagnose, treat, cure, or prevent any disease. If you have any medical conditions, consult a healthcare professional before use.

HOW TO USE:

Kit Includes 6 bottles of UTRUSTME Micro infusion Serum and 1 Derma Roller Head.**Treatment Time:** 5 minutes

Recommended Use: 1 treatment per week for 6 weeks.

How to Micro-Infuse:

- Cleanse your face, remove all makeup and ensure there is no dirt, oil or stray hairs on your face that could increase the risk of infection or irritation.
- Wash your hands and ensure you have a sterile environment for the treatment.
- Remove the cap from a bottle of micro-infusion serum
- Place the micro-needle head onto the serum bottle. **NOTE: Spray or bathe the micro-needle head with alcohol to sanitise it after and before each recurring use. Failure to do this may cause infection.**
- Hold the bottle upside down for 30-60 seconds to allow the serum to drip onto the micro-needle head. If it's easier, shake the bottle. You may need to repeat this a few times in the treatment to ensure the needles are covered in serum.
- Lightly roll the micro-needle head over the skin in an upward direction (following the diagram on the next page) to cover the forehead, cheeks, chin and jawline. Avoid using it on the neck, as the skin is very fragile.
- Roll 2-3 times over each section until you have covered your face, ensuring to keep the needles wet with serum as you go. This treatment should take no longer than 5 minutes.
- After the treatment, you will most likely have serum remaining in the bottle. Remove the micro-needle head and apply the remaining serum to your face by hand.
- Repeat weekly for 6 weeks.

Which direction to roll the micro-needle head:



INGREDIENTS & LID REMOVAL:

Our micro infusion serum has been expertly formulated to give you the ultimate treatment; it contains the following ingredients:

Pure water, Hyaluronic acid, Hansheng gum, Allantoin, Biological sugar gum, BSASM, Oligopeptide-1, P-hydroxy acetophenone, 1,2-hexanediol, SM Nicotinamide (VB3), HP, MOSK 100, Hexapeptide-8, Conopeptide, Bosein (powder), Ikedoin, Ergothione.



SAFETY GUIDELINES:

Things TO DO:

- Seek advice from a dermatologist if you have any skin conditions prior to using the product.
- Ensure your hands are thoroughly cleansed before and during the treatment process. If you prefer an added layer of protection, you can opt for gloves.
- It's highly recommended to conduct a patch test on a discreet area of your skin before treatment. This test is crucial as it helps you anticipate potential reactions and recovery time, ensuring a safe and effective treatment process.
- Prior to use, sanitise both your hands and the device.
- In the event of any bleeding, cleanse your face with water.
- Apply a high-SPF sunscreen post-treatment.
- Store the device in the provided box for optimal hygiene and safety. After each treatment, sterilize the device thoroughly with isopropyl alcohol or a similar disinfecting agent.

Things to AVOID:

- Improper sterilisation of the needle head can lead to infection, irritation, or cross-contamination. Always disinfect before and after use with isopropyl alcohol.
- Avoid applying topical retinoids for at least a day before undergoing treatment.
- Limit sun exposure, especially in the three days leading up to your treatment.
- When treating areas with active blemishes, approach as closely as possible but refrain from going over any active breakouts.
- Do not proceed with treatment if your skin is already irritated or inflamed.
- Do not use if the needle head is damaged or compromised. Broken needles may cause skin trauma or infection.
- Refrain from repeatedly treating the same area of your face.
- Apply gentle pressure and avoid using excessive force during treatment.
- Avoid using comedogenic serums or oils during or after treatment at home.
- Exercise caution when treating raised scars such as keloids.
- Refrain from using active skincare products until at least three days after treatment.
- Minimize sun exposure to the treated area for up to 72 hours following the treatment session.

POST TREATMENT

Things TO DO:

- **CLEANSE AND HYDRATE** The day following your treatment, cleanse your face with a mild cleanser and apply a hydrating moisturizer to replenish your skin's moisture.
- **MINIMISE MAKEUP USAGE** It's advisable to refrain from using makeup for as long as possible after treatment, ideally for a minimum of 24 hours. If using makeup cannot be avoided, opt for products designed for post-treatment care or formulated for sensitive skin to reduce the risk of irritation.
- **AVOID HARSH INGREDIENTS** After using the UTRUSTME Micro Infusion Kit, be mindful of the products you apply to your skin, as they will absorb more readily. Opt for gentle formulations without harsh ingredients, specifically avoiding acids to prevent any potential discomfort.
- **STEER CLEAR OF WATER ACTIVITIES** To protect your skin's integrity, it's best to avoid swimming, saunas, and steam rooms for at least 24 hours following your treatment.
- **MINIMISE REDNESS** Following the procedure, some users may experience minor redness and irritation, which should resolve quickly. You may also notice mild swelling, peeling, or flaking. Keeping your skin adequately moisturised can help reduce shedding and alleviate any feelings of tightness.
- **LIMIT SWEATING** For the initial 24 hours post-treatment, refrain from engaging in vigorous exercise or activities that induce sweating. Gym environments, in particular, harbour bacteria and can potentially trigger adverse reactions.
- **HANDS OFF** Avoid touching your face, especially with unwashed or unsanitised hands, to minimise the risk of introducing bacteria or causing irritation.

POST TREATMENT

Things to AVOID:

- **Refrain** from washing your face immediately following the treatment.
- **Avoid** using makeup for at least 24 hours after the procedure. If makeup application is necessary, opt for products designed for sensitive skin or post-treatment use to minimise the risk of irritation.
- **Minimise** direct sun exposure for 24-48 hours post-treatment. Apply sunscreen with high SPF before stepping outdoors and limit prolonged sun exposure. Choose a non-comedogenic sunscreen for sun protection.
- **To prevent** potential irritation, refrain from applying skin care products containing fragrances, active ingredients, or acids such as Alpha-Hydroxy Acids, Beta-Hydroxy Acids, Retinol (Vitamin A), Vitamin C, or other potent skincare products like scrubs or toners for at least 2-3 days following the treatment.
- **Avoid** using exfoliating products for 2-3 days after the treatment as they may exacerbate skin sensitivity during this period.

WARNINGS:

- **Keep** out of reach of children; some small parts may pose a choking hazard.
- **Discontinue** use if any concerns arise and speak to your doctor. Only use and apply products after treatment that are not compromised or contaminated.
- **A patch test** is required before proceeding with the full treatment.
- **If you** are unsure if the product suits you, please seek advice from your doctor.
- **Take care** when using the device if you have a history of post-inflammatory hyperpigmentation. People with a history of this are more likely to pigment after microneedling.
- **Before use**, ensure your skin is thoroughly cleaned and free from products or contaminants.
- **Do not** perform more than two treatments per week.
- **Ensure** that the needle heads do NOT come into contact with contaminated or unclean objects such as clothing or benchtops
- **Avoid** UV radiation and sunlight after treatment.
- **This product** requires the use of needle heads, which, under certain circumstances, pose the risk of a sharps injury. Injuries from contaminated needles in needleheads may lead to the transmission of diseases.
- **Take care** when handling needle heads, as there is a risk of a sharps injury. If all needles are NOT fully intact, the needle head is faulty and must be immediately replaced.
- **This product** contains a safety protection feature built into the derma roller head to prevent the needle from going any deeper than 0.5 mm.
- **Do NOT** attempt to reuse the needle head without soaking or spraying with isopropyl alcohol to sanitise the head first.
- **Do NOT** touch the needles and tip at any time. Contact with the needles may lead to a sharps injury or contamination.
- **Use ONLY** needle heads on the device that has been provided.
- **This** treatment is not intended to diagnose, treat, cure, or prevent any disease. If you have any skin conditions or are pregnant or breastfeeding, please consult with your doctor prior to treatment.
- **If you** have metal allergies or if you have any medical concerns, please consult your medical professional before using this treatment device. Any harmful consequence resulting from misuse or application to other body areas, not following instructions or any other improper applications is neither the responsibility of UTRUSTME or its affiliates.

CONTACT

Sensitivity Assessment:

Prior to commencing your initial treatment, it is essential to perform a sensitivity assessment.

1. Remove the cap off a serum bottle and place a small amount on your inner arm.
2. Allow a minimum of 3 hours to elapse. If your inner arm exhibits irritation or discomfort, this indicates that the serum may not be suitable for you.

Contact:

If you have any further questions or doubts about using the **UTRUSTME** Micro Infusion Kit, please get in touch with us using the details below, and our support team will assist you.

If you have any medical concerns or are unsure if the device suits you, please consult a medical professional before use.

If you experience any unexpected side effects, skin damage, or infections, discontinue use and report it to our support team.

If serious reactions occur, seek medical attention immediately.