

Disposable Medical Face Mask

Epidemic emergency products

Non-sterile

【Product name】

Disposable medical face mask

【Model, specifications】

Model: flat ear-hook type; Specifications: 17.5cm * 9.5cm.

【Structure and composition】

The disposable medical face masks are mainly composed of a mask body, a mask, and a nose clip. The mask body is made of non-woven fabric (inner and outer layer) and melt-blown cloth (middle layer).

【Application】

It is suitable for covering the user's mouth, nose and jaw. It is used in ordinary medical environments to wear and block the mouth and nasal cavity from exhaling or spraying pollutants.

【Product performance】

The main performance of the product complies with the CYT / T 0969-2013 Disposable Medical Masks standard.

【Product wearing method】

1. In general, the colored side/the side with folds facing downwards of the face mask should face outwards with the metallic strip uppermost
2. Mould the metallic strip over nose bridge and the face mask should fit snugly over the face
3. Extend the face mask to fully cover mouth, nose and chin

【Contraindications, precautions, warnings and tips】

- 1) This product cannot produce oxygen and cannot be used in an anoxic environment;
- 2) Do not use this product when facial beard, hair, or other conditions can affect the adhesion between the face and the mask;
- 3) This product is only for single use by a single patient, and repeated use is prohibited;
- 4) Check whether the packaging is intact before use, and the single package is damaged and prohibited to use
- 5) Please note that pulling the mask band too strongly will cause the mask band to break;
- 6) If your mouth, nose, ears, or the surface of the skin that comes into contact with this product when you expect to use this product are injured, treat it before using this product;
- 7) After using this product, it should be disposed of in accordance with hospital regulations or environmental protection department requirements;

【Storage and transportation conditions】

This product should be stored in a dry, dark, heat-proof, ventilated, and non-corrosive gas place; avoid direct sunlight and rain and snow during transportation.

【Production date / batch number】

Printed on the package

【Validity Period】

2 years

【Interpretation of graphics, symbols, abbreviations, etc. used in medical device labels】

1. For 1 time use only.
2. Please read the manual carefully.
3. Keep away from rain
4. Keep out of the direct sun

【Implementation standards】 YY/T0969-2013 （disposable medical face mask）