## MOIST & HEAT PAD

Ergonomic Digital Hand Control with an LED indicator allows operator to program, and monitor the temperature and treatment time.



Far-Infrared

its own moisture

OO-053 NO WATER NEEDED 20"X20" For the automatically makes joint of shoulder.

relieve pain caused by muscle spasm or stiffness.

- No water needed, automatically makes its own moisture.
- Temporary effectively relieve of arthritic and rheumatic with the use of moist pad.
- Relieve the injury and pain after sporting or working caused.



▲ OO-052 14"X27" for large area of body, back, entire spine. ▲ OO-054 14"X20" for large area of body, back, entire spine.







OO-055 14"X14" For knee, shoulder, small area of body.





00-064 4"X20" Special size for neck or hollow part of body.



Hand Control with an LED indicator.



Moisture is drawn from humidity in the air to produce steam.

## NO WATER NEEDED

automatically makes its own moisture



▲ OO-046 14"X27" for large area of body, back, entire spine.



14"X14" For knee , shoulder , small area of body .



4"X20" Special size for neck or hollow part of body.



- Moist heat is especially useful in relieve pain caused by muscle spasm or stiffness.
- No water needed, automatically makes its own moisture.
- Temporary effectively relieve of arthritic and rheumatic with the use of moist pad.
- Relieve the injury and pain after sporting or working caused.



Switch with rolling adjustment.



Combination of ceramic far-infrared elements.



Moisture is drawn from humidity in the air to produce steam.





Code	Size	Voltage		Power Consumption
00-056	30 cm x 36 cm	Ac 110	60HZ	55 W
00-058	30 cm x 60 cm	Ac 110	60HZ	55 W
00-059	30 cm x 36 cm	Ac 220	50HZ	55 W
00-060	30 cm x 60 cm	Ac 220	50HZ	55 W

## **GOOD** assistant for home care

- Keeps you warm.
- Eases your pain.
- Relaxes your body.
- Helps your body's fluids circulation.



Heating application can ease and cure arthritis as well as rheumatism.



Heating application provides excellent treatment.



As soon as you feel pain in your muscles and joints, use the Dry & Moist Heating Pad to ease the hurt.



Heating application can effectively cure back pain.

