

Medical prevention product brochur

XIAMEN ECO PACK Co.,Ltd





XIAMEN ECO PACK CO LTD

Product introduction 1



EC014

PRODUCT NAME: Daily protective mask

COMPONENTS : This product consists of a mask body, a nose clip, and a mask band

SCOPE AND USE:

Suitable for filtering particulate matter in air pollution environment in daily life

INSTRUCTION: Wear on the nose and mouth

WARNING:

- 1.Period of validity is 3 years, please use it within this limited period.
- 2.For one-time use only
- 3.Please do not use it when package in damaged condition
- 4.Using with caution for those who are allergic to nonwovens
- 5.Please refer to the instruction before using.

CONTRAINDICATION: none

SPECIAL STORAGE CONDITIONS AND METHODS:

Store in a room which is in a humidity not exceeding 80%, non-corrosive gas and well ventilation, avoid high temperature.

VALIDITY: validity is 3 years under specified conditions storage

IMPLEMENTATION STANDARDS: GB/T 32610



XIAMEN ECO PACK CO LTD

Product introduction 2



EC002

PRODUCT NAME: Daily protective mask

COMPONENTS : This product consists of a mask body, a nose clip, and a mask band

SCOPE AND USE :

Suitable for filtering particulate matter in air pollution environment in daily life

INSTRUCTION: Wear on the nose and mouth

WARNING:

- 1.Period of validity is 3 years, please use it within this limited period.
- 2.For one-time use only
- 3.Please do not use it when package in damaged condition
- 4.Using with caution for those who are allergic to nonwovens
- 5.Please refer to the instruction before using.

CONTRAINDICATION: none

SPECIAL STORAGE CONDITIONS AND METHODS:

Store in a room which is in a humidity not exceeding 80%, non-corrosive gas and well ventilation, avoid high temperature.

VALIDITY: validity is 3 years under specified conditions storage

IMPLEMENTATION STANDARDS: GB/T 32610-



XIAMEN ECO PACK CO LTD

Product introduction 3



EC013

PRODUCT NAME: N95-Protective mask

COMPONENTS : This product consists of a mask body, a nose clip, and a mask band

SCOPE AND USE:

Suitable for filtering particulate matter in air pollution environment in daily life

INSTRUCTION: Wear on the nose and mouth

WARNING:

- 1.Period of validity is 3 years, please use it within this limited period.
- 2.For one-time use only
- 3.Please do not use it when package in damaged condition
- 4.Using with caution for those who are allergic to nonwovens
- 5.Please refer to the instruction before using.

CONTRAINDICATION: none

SPECIAL STORAGE CONDITIONS AND METHODS:

Store in a room which is in a humidity not exceeding 80%, non-corrosive gas and well ventilation, avoid high temperature.

VALIDITY: validity is 3 years under specified conditions storage

IMPLEMENTATION STANDARDS: GB 2626-2006



XIAMEN ECO PACK CO LTD

Package 1



x 50



x 40



Color box (50 Pieces)
Material is 350g gray paper
1. Size:195*105*80mm,
Tolerance:±1mm
2. 4C: Printing

Outer case (40 Color boxes)
The material is K=K
1. Size: 545*410*345mm,
Tolerance: ±5mm
2. Bursting strength: 26.5kg/cm
Bearing strength : 0.05kgf/cm

The above is for reference only, subject to final delivery.



XIAMEN ECO PACK CO LTD

Package 3



x 50



x 40



Color box (50 Pieces)
Material is 350g gray paper
1. Size: 140*100*80mm,
Tolerance: ± 1 mm
2. 4C: Printing

Outer case (40 Color boxes)
The material is K=K
1. Size: 425*300*430mm
Tolerance: ± 5 mm
2. Bursting strength: 26.5kg/cm
Bearing strength : 0.05kgf/cm

The above is for reference only, subject to final delivery.



EOP19

x 10
→



Color box (10 Pieces)
Material is 350g gray paper
1. Size:108*157mm,
Tolerance:±1mm
2. 4C: Printing

x 400
→



Outer case (400 Color boxes)
The material is K=K
1. Size: 370*370*550mm
Tolerance: ±5mm
2. Bursting strength: 26.5kg/cm
Bearing strength : 0.05kgf/cm

The above is for reference only, subject to final delivery.