

Foot Massager 002

[Product Function Hotspots]

1. Large area vibration moxibustion on the sole of the foot, from the finger to the heel.
2. The vibration intensity is infinitely variable and can be adjusted at will.
3. The particles protruding from the surface of the vibrating part stimulate the local community and activate the meridians.
4. When using, vibration, magnetic therapy and hyperthermia can be combined into one or used separately, and can be selected at will.

[Main technical indicators]

Rated voltage: AC 220V, frequency 50Hz

Total input power 3W

Amplitude $\pm 2.0\text{mm}$

Thermal temperature above 40 °C

[Precautions]

1. It is forbidden for pregnant women and children to use it;
2. Suspend use within 30 minutes after meals
3. Use twice every morning and evening for 15 minutes each time
4. Cut off the power after use
5. When unplugging the power supply from the socket, please hold the plug and pull it out. Do not pull the power cord directly
6. Please use a soft cloth to clean the volatile chemicals of the warrior figurines.
7. Prevent water from entering the machine to avoid failure.
8. The vibration power is 16-22W, the rated working time is 15 minutes, the total input power is 3W, and the thermal temperature is above 40 °C.
9. It is recommended that each use should not exceed 15 minutes. When using again, please keep an interval of more than 15 minutes

Packing quantity: 16 sets

Outer box size: 68 * 37 * 67cm

Gross weight of one set: 1.5kg

Gross weight of one box: 25kg

Volume weight of one box: 34kg

Color box size: 32.5 * 36 * 9cm

Product size: 34 * 31 * 5cm