


Reliable Alarm Activation at Any Corner



Raise an alarm by pulling the attached cord or pressing the button



IPX4 Waterproof standard



Wall-mount design



PCU-8

Pull Cord Unit

The Pull Cord Unit PCU-8 is used to activate an alarm when the cord is pulled or when the button is pressed, allowing users to quickly activate an alarm for emergency help. The added button feature enhances usability, providing an alternative method for triggering emergency assistance.

The wall-mounted PCU-8 is a crucial safety device, particularly useful when the user is not wearing a personal trigger. Once triggered, the alert is transmitted via the connected medical alarm to monitoring centers or care personnel. Rated at IPX4, it is resilient against water splashes, ensuring reliability in various indoor environments.

PCU-8

Pull Cord Unit



Other Features

- Raise an alarm by pulling the attached cord or pressing the button
- IPX4 Waterproof standard
- Extensive RF communication range & faster signal transmission
- Wall-mount design
- Low battery detection
- Supervision signals to check system integrity

Specifications

PCU-8 Pull Cord Unit

Frequency	919 MHz
Power Source	1.5V, Alkaline batteries x2
Battery Life	10 years*
Waterproof Level	IPX4
Cord Length	91±2cm
Operating Temperature	-10°C to 45°C (14°F to 113°F)
Operating Humidity	Up to 85% non-condensing
Dimensions	99mm x 85mm x 29mm

**Note: Actual battery life varies with product settings, conditions, operating environment, and usage.*