

sleep aid HSA009

Product principle of ear clip type sleep aid: mainly through CES therapy, also known as "transcranial micro-current stimulation therapy", the micro-current stimulates central nerve cells through the brain, reduces anxiety and increases attention, and changes the sensory signal to make the brain ALPHA (alpha ray) wave rise, making people feel relaxed, regulating mood, promoting sleep, and improving sleep quality.

Features: 1: clip ear appearance, comfortable ear clip, no feeling, no falling

2: Three modes, comfort mode and stimulation mode. Strengthen mode selection.

3: Extra long endurance. 4 - 5 hours at a time, available for 2 weeks (14 days)

Input voltage: DC 5V/1A

Rated power: 1.5W

Working current: 0.2A-0.5A

Function mode: θ Theta wave/ α Alpha wave

Pulse gear: 9-gear adjustment

Charging time: 5H

Duration of use: 20-60 times

Battery capacity: 1500mAh/3.7V

Product size: 136*65*18mm

Color box size: 187*146*72mm

Weight: 0.2KG