

Mini CES Sleep Aid 1115

Host dimensions: 65*30*23mm

Material: ABS resin + PU leather + conductive silicone

Operating voltage/current: 3.3V/50mA

Working time: 15 minutes per session

Built-in lithium battery: 80mAh

Operating frequency: 1-100Hz

Smart chips enable more precise control

Physical microwave pulses, non-invasive, painless, and worry-free

Massage for sleep aid—apply where sore

Alleviate insomnia, improve sleep quality, and promote sleep scientifically

10 Massage Modes: Acupuncture and Moxibustion, Tapping, Kneading, Pressing, Massage, Percussion, Beating, Rolling, Stretching, Finger Pressure

20-level intensity adjustment

1-8 levels: Low intensity: Gentle current for sensitive areas

Levels 9-14: Moderate intensity, suitable for most individuals with mild to moderate pain

15-20 levels: Strong intensity, suitable for most patients with moderate to severe pain

Contraindications:

1. Contraindicated for acute patients with cerebral hemorrhage
2. Contraindicated for patients with epilepsy or pregnant women
3. Victims who have suffered a major blow or are just recovering from it should avoid it
4. Contraindicated in patients with severe depression and recent stroke
5. Contraindicated for individuals using cardiac pacemakers or other electrotherapy devices
6. Patients with earlobe damage or non-healed earlobe wounds are prohibited from use
7. If discomfort occurs during use, discontinue immediately

Weight: 0.35KG