

Cervical vertebra massage instrument M008

Product material: ABS

Rated voltage: 5V

Rated power: 12W

Color box size: 16 * 14.5 * 5.2cm

Weight: 0.18KG

Outer box size (60PCS): 52 * 31 * 50CM

Weight: 10.8KG (volume weight: 16.2KG)

The trinity of low-frequency electric pulse, magnetic effect and circular traction forms an efficient composite energy field,

TEMS technology: use the pulse electric frequency to transmit to the deep skin, simulate a variety of massage, and bring a three-layer deep massage experience

Gear 12: gear 1-3 is gentle to relieve muscle stiffness, gear 4-6 is medium to relieve lactic acid accumulation, gear 6-9 is strong to relieve muscle soreness, gear 9-12 is strong to improve muscle soreness and discomfort.

Constant temperature hot compress: one button heating, built-in new temperature controller, to achieve 3 seconds of rapid heating, further eliminate fatigue.

Ergonomic design: It adopts 3D intelligent fitting technology, which can easily adjust the electrode according to the curve of your neck.