


<p><b>KCA435</b></p>	<p>Ultrasoun d Cavitation Body Slimming Massager</p>	
----------------------	--	---

Ultrasonic Frequency:100K times/s  
Unit size: 103\*103\*111mm  
Operating Voltage:8.4V  
Operating Current:0.6-1A  
Power:8.4W  
levels:3 levels  
Charging time: 2 hours  
Material: ABS+Zinc alloy  
N.W.:326g  
Technology: RF cavitation, EMS,  
Vibration , red/ blue leds ,  
Function: Accelerate Metabolism,  
Enhance Tissue Regeneration  
Process, Dilute Dullness, Anti Aging,  
Lose Weight, Slim Body

Rejuvenating Mode:  
RF+infrared+red light  
Fat Bruning Mode: Ultrasound  
+vibration+infrared+red light  
Fat Breaking Mode:RF+ Ultrasound  
+vibration+infrared+red light  
Shaping  
Mode: RF+ EMS+Ultrasound  
+infrared+blue light

Ultrasound: The fat-burning machine can remove fat sacs deep in the skin, reduce stubborn fat, achieve the effect of weight loss and body shaping, and shape the abdomen, arms, waist, legs, and buttocks. This is the most correct decision for you to choose a body sculpting massager.

Infrared+Red Led light therapy stimulates cellular activity to help minimize fine lines and wrinkles

EMS therapy sends gentle waves through the skin, tissues, and muscles to stimulate collagen and elastin to help firm and tighten

Bi-Polar RF stimulates collagen and penetrates deep into the skin tissue to help the body eliminate fat cells

Vibration: The vibration of the fat-burning machine shrinks collagen fibers, accelerates the rate of fat decomposition, and effectively reduces stubborn fat.